




Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 1 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

1
2
3
4
5
6
7
8
9
10


**LabCyte CORNEA-MODEL24
EYE IRRITATION TEST
OPERATION PROTOCOL
Ver. 2.7**

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 2 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

11	TABLE OF CONTENTS	
12		
13	1. RATIONALE AND BACKGROUND	4
14	1.1 The LabCyte CORNEA-MODEL24 EYE IRRITATION TEST	4
15	1.2 BACKGROUND OF LabCyte24 EIT	4
16	1.3 BASIS OF THE TEST METHOD.....	5
17	1.4 LIMITATION OF THE METHOD.....	5
18	1.5 BRIEF BASIC PROCEDURE.....	5
19	1.5.1 LabCyte CORNEA-MODEL24	6
20	1.5.1.1 QUALITY CONTROL OF THE TEST SYSTEMS	6
21	1.5.1.2 PRECAUTIONS.....	7
22	1.5.2 ASSAY QUALITY CONTROL	7
23	1.5.2.1 ASSAY ACCEPTANCE CRITERION 1: NEGATIVE CONTROL	7
24	1.5.2.2 ASSAY ACCEPTANCE CRITERION 2: POSITIVE CONTROL.....	7
25	1.5.2.3 ASSAY ACCEPTANCE CRITERION 3: STANDARD DEVIATION (SD)	8
26	1.6 DATA INTERPRETATION PROCEDURE (PREDICTION MODEL).....	8
27	2. MATERIALS.....	9
28	2.1 LabCyte CORNEA-MODEL24.....	9
29	2.1.1 LabCyte CORNEA-MODEL24 KIT COMPONENTS	9
30	2.1.2 SHIPMENT OF LabCyte CORNEA-MODEL24	9
31	2.1.3 INSTRUCTIONS FOR USE OF LabCyte CORNEA-MODEL24.....	10
32	2.2 CONSUMABLES.....	10
33	2.3 OTHERS	11
34	2.3.1 EQUIPMENT/INSTRUMENTS	11
35	3. TEST METHOD.....	12
36	3.1 PREPARATIONS.....	12
37	3.1.1 POSITIVE CONTROL.....	12
38	3.1.2 NEGATIVE CONTROL.....	12
39	3.1.3 POLY WASH BOTTLE FOR PBS	12
40	3.2 TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT	13
41	3.2.1 DETECTION OF THE CHEMICALS THAT STAIN THE TISSUE	13
42	3.2.1.1 STEP 1 (PRELIMINARY TEST).....	13
43	3.2.1.2 STEP 2 (FUNCTIONAL CHECK ON VIABLE TISSUE).....	13
44	3.2.2 DETECTION OF CHEMICALS THAT DIRECTLY REDUCE WST-8	14
45	3.2.2.1 STEP 3 (PRELIMINARY TEST).....	14
46	3.2.2.2 STEP 2 (FUNCTIONAL CHECK ON FREEZE-KILLED TISSUE)	15
47	3.3 EXECUTION OF THE TEST	17
48	3.3.1 PREPARATION OF LabCyte CORNEA-MODEL24 (DAY -1).....	17
49	3.3.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1).....	18
50	3.3.2.1 PREPARATION OF WELLS FOR POST-INCUBATION (3 RD ROW)	18
51	3.3.2.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING.....	19
52	3.3.2.3 POST-EXPOSURE INCUBATION	21
53	3.3.3 APPLICATION OF SOLID TEST CHEMICALS AND RINSING (DAY 0~1).....	22
54	3.3.3.1 PREPARATION OF WELLS FOR SOLID TEST CHEMICAL APPLICATION (2 ND ROW) AND FOR	
55	POST-INCUBATION (3 RD ROW).....	22
56	3.3.3.2 APPLICATION OF SOLID TEST CHEMICALS.....	22

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 3 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

57	3.3.3.3	RINSING OF SOLID TEST CHEMICALS.....	23
58	3.3.4	WST-8 ASSAY (DAY 1).....	25
59	3.3.4.1	PREPARATION OF WELLS FOR WST-8 ASSAY.....	25
60	3.3.4.2	WST-8 REACTION.....	25
61	3.3.4.3	SAMPLING THE REACTED WST-8 MEDIUM.....	26
62	3.3.4.4	OPTICAL DENSITY MEASUREMENTS OF THE REACTION MEDIUM.....	27
63	4.	ASSESSMENT.....	29
64	4.1	CONDITIONS FOR A VALID TEST.....	29
65	4.2	ASSAY CRITERIA.....	29
66	5.	References.....	30
67	MDS 1-1:	RECEIPT OF LabCyte CORNEA-MODEL24 (2.1.2).....	31
68	MDS 1-2:	TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8	
69		ENDPOINT STEP-1 (3.2.1.1).....	32
70	MDS 1-3:	TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8	
71		ENDPOINT STEP-1 (3.2.2.1).....	33
72	MDS 1-4:	PREPARATION OF FREEZE KILLED TISSUE (3.2.2).....	34
73	MDS 2:	PREPARATION OF LabCyte CORNEA-MODEL24 (3.3.1).....	35
74	MDS 3-1(LIQUID):	APPLICATION OF LIQUID TEST CHEMICALS, RINSING AND	
75		POST-INCUBATION (3.3.2).....	36
76	MDS 3-2(SOLID):	APPLICATION OF SOLID TEST CHEMICALS AND RINSING (3.3.3)	
77		38
78	MDS 4-1:	WST-8 ASSAY (3.3.4.1, 3.3.4.2).....	40
79	MDS 4-2:	SAMPLING OF REACTING WST-8 SOLUTION (3.3.4.3) AND	
80		MEASUREMENT (3.3.4.4).....	41
81		REVISION HISTORY.....	42
82			

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 4 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

83 **1. RATIONALE AND BACKGROUND**

84 **1.1 The LabCyte CORNEA-MODEL24 EYE IRRITATION TEST**

85 The LabCyte CORNEA-MODEL24 eye irritation test (LabCyte24 EIT) is designed to identify test
86 chemicals that cause acute eye irritation by measuring cytotoxic effects using the
87 2-(2-methoxy-4-nitrophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt
88 (WST-8) assay on the reconstructed human corneal epithelial (RhCE) model. The LabCyte24 EIT is not a
89 kit, but LabCyte CORNEA-MODEL24 tissues are commercially available at a minimum of 24 LabCyte
90 CORNEA-MODEL24 tissues per order.


91

92 **1.2 BACKGROUND OF LabCyte24 EIT**

93 Assessment of ocular irritation is an essential part of early testing procedures for the evaluation and
94 hazard classification of substances. Therefore, it plays an important role in the safety evaluation of
95 general consumer products and materials. Novel substances used in consumer products must undergo
96 comprehensive toxicological evaluation for eye irritation and a variety of other adverse outcomes. To
97 date, the *in vivo* Draize eye test has been the international standard assay for acute ocular toxicity
98 evaluation (irritation and corrosion) and is described (including optimizations and refinements) in OECD
99 Test Guideline 405. However, the use of this test has been questioned and strongly criticized for ethical
100 concerns related to animal welfare, because it is painful to the rabbits. Thus, alternative strategies and
101 tests are urgently required in order to evaluate the eye irritation potential of new chemicals. Corneal
102 epithelial cells on the surface of the eye are to first to be exposed to substances and have been widely
103 studied for links to the biological role of tissue and gene regulation. Three-dimensional RhCE models are
104 useful as a multilayered, standardized tissue that mimics the human corneal epithelium.

105 The LabCyte24 EIT was developed as a replacement for the Draize eye irritation test. The Draize
106 scoring system is heavily weighted towards corneal damage (80 out of 110 total score), because
107 irreversible damage to the cornea can lead to blindness. Since damage to the cornea is so crucial for
108 human health, corneal tissue can be considered a useful tool for the development of *in vitro* eye irritation
109 testing.

110

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 5 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

111 **1.3 BASIS OF THE TEST METHOD**

112 Damage induced by eye irritants generally progresses from the corneal epithelium through the stroma
113 and potentially to the endothelium, and the LabCyte24 EIT is able to provide information on the first
114 stages of this progression. Irritants damage cells while penetrating the corneal epithelium layer, and the
115 cytotoxic progression can be estimated by analyzing cell viability in the LabCyte CORNEA-MODEL24
116 tissue using standardized methods. Although the tissues represent only the corneal epithelium, (very mild
117 responses would also be reflective of some conjunctival irritation), it can be used to estimate deeper
118 damage as far as the stroma by the analyzing cell viability.

119 The relative viability of the tissue exposed to a test chemical is measured using the WST-8 assay
120 immediately after exposure and again after a post-exposure period. A viability of 40% of the negative
121 control value was used as the cutoff in identifying test chemicals as an irritant (GHS Category 1 or 2) or
122 an non-irritant (GHS No Category). Some culture environments might allow the detection of very small
123 quantities of cytokines secreted by the corneal epithelial tissue in response to topical application of test
124 chemicals.

125


126 **1.4 LIMITATION OF THE METHOD**

127 One limitation of this assay method is potential interference of the test chemical with the WST-8
128 endpoint. A colored test chemical or one that directly reduces WST-8 (and thereby mimics dehydrogenase
129 activity of the cellular mitochondria) could interfere with the WST-8 endpoint. This is only a problem,
130 however, if there are significant residual levels of the test chemical on or in the tissue at the time of the
131 WST-8. Although this is considered unlikely, if it did happen, both the actual metabolic WST-8 reduction
132 and the quasi direct WST-8 reduction by a colored test chemical can be quantified using the procedure
133 described in **Section 3.2 “TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH**
134 **WST-8 ENDPOINT”**.

135

136 **1.5 BRIEF BASIC PROCEDURE**

137 On the day of receipt, LabCyte CORNEA-MODEL24 tissues are incubated overnight to release
138 transport–stress-related compounds and debris.

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 6 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

139 For liquid test chemicals, tissues are topically exposed for 1 minute. Preferably, three tissues are used
140 for each test chemical as well as for the positive and negative controls. After exposure, tissues are
141 thoroughly rinsed and blotted to remove the test chemical, then transferred to fresh medium and
142 post-incubated for 24 hours. For solid test chemicals, tissues are exposed for 24 hours but are not
143 post-incubated.


144 After post-incubation of tissue exposed to liquid test chemicals or after exposure of tissue exposed to
145 solid test chemicals, the tissues are each transferred to a well containing the WST-8 medium in a 1:10
146 dilution with Earle balanced salt solution (EBSS). After a four-hour WST-8 incubation, the orange
147 water-soluble formazan salt is formed in the WST-8 medium by cellular mitochondria and optical density
148 (OD) of the WST-8 medium is measured using a spectrophotometer at 450 nm and 650 nm as reference.
149 Relative cell viability is calculated for each tissue as % of the mean of the negative control tissues. Test
150 chemicals that produce a relative cell viability below 40% of the negative control are predicted to be
151 irritants.

153 **1.5.1 LabCyte CORNEA-MODEL24**

154 The LabCyte CORNEA-MODEL24 is a commercially available RhCE model produced by Japan
155 Tissue Engineering Co. Ltd. It comprises normal human corneal epithelial cells that are derived from a
156 human cornea. The cells are cultured with 3T3-J2 cells as a feeder layer in order to expand them while
157 maintaining their phenotype (Rheinwald and Green, 1975; Green, 1978). Reconstruction of the human
158 cultured corneal epithelial tissue is achieved by cultivating and proliferating the corneal epithelial cells
159 on an inert filter substrate with a surface area of 0.3 cm² at an air-liquid interface for 13 days using an
160 optimized medium containing 5% fetal bovine serum. This results in the formation of a multilayer
161 structure comprising a fully differentiated corneal epithelium with features mimicking those of a normal
162 human corneal epithelium. For delivery, LabCyte CORNEA-MODEL24 tissues are embedded in an
163 agarose gel containing a nutrient solution and shipped in 24-well plates.

165 **1.5.1.1 QUALITY CONTROL OF THE TEST SYSTEMS**

166 LabCyte CORNEA-MODEL24 tissue is manufactured in accordance with well-defined quality

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 7 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

167 assurance procedures. Each production batch was comes with quality control documentation that
168 identifies storage conditions, RhCE instructions for use, lot number and origin, histology (demonstration
169 of human multilayered corneal epithelial-like structure), cell viability, and barrier function integrity
170 ($0.1 \leq IC_{50} \leq 0.4$).

171
172

173 **1.5.1.2 PRECAUTIONS**

174 Corneal epithelial cells are taken from healthy donors who are free of HIV or hepatitis. Nevertheless,
175 always adhere to the following procedures for the handling of biological materials:

- 176 (a) Always wear gloves during handling of the eye and kit components.
- 177 (b) All corneal epithelial tissue as well as all material and media that came in contact with it should be
178 decontaminated prior to disposal using special containers or autoclaving.

179

180 **1.5.2 ASSAY QUALITY CONTROL**

181

182 **1.5.2.1 ASSAY ACCEPTANCE CRITERION 1: NEGATIVE CONTROL**

183 The absolute OD of the negative control for either liquid or solid test chemicals (NC-L or NC-S)
184 tissues (treated with sterile PBS for liquid test chemicals or untreated for solid test chemicals) in the
185 WST-8 assay is an indicator of tissue viability obtained in the testing laboratory after shipping and
186 storing procedures and under specific conditions of use.

187 **$0.5 \leq \text{Mean OD (A450/650) measured value} \leq 1.6$**


188

189 **1.5.2.2 ASSAY ACCEPTANCE CRITERION 2: POSITIVE CONTROL**

190 Ethanol is used as the positive control (PC) for liquid test chemicals (PC-L) and is tested concurrently
191 with the liquid test chemicals. Lauric acid is used as the PC for solid test chemicals (PC-S) and is tested
192 concurrently with the solid test chemicals.

193 Concurrent here means that the PC-L and the PC-S are to be tested for each run.

194 **Mean tissue viability $\leq 40\%$**

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 8 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

195

196 **1.5.2.3 ASSAY ACCEPTANCE CRITERION 3: STANDARD DEVIATION (SD)**

197 Since eye irritation potential is predicted from the mean viability of three individual tissues, the
198 variability of tissue replicates must kept at an acceptably low level.

199 **Standard Deviation (SD) of tissue viability of three identically treated replicates for the negative**
200 **control, positive control, and test chemicals \leq 18%**

201

202 **1.6 DATA INTERPRETATION PROCEDURE (PREDICTION MODEL)**

203 The United Nations Globally Harmonized System (GHS) (United Nations, 2003) was used as a
204 reference for the in vivo eye irritation classification of test chemicals.

205 For the purpose of this EIT, an eye irritant is defined as a substance that induces reversible ocular
206 lesions after administration to rabbits.


207 According to GHS classifications, a substance is an irritant (Category 1 or 2) if the mean relative
208 viability of three individual tissues exposed to the test chemical is falls below 40% of the mean viability
209 of the negative control. **(Refer to Table 1.)**

210 **Table 1** Prediction model of LabCyte24 EIT

<i>In vitro</i> results	Prediction
Tissue viability is \leq 40%	Category 1 or 2
Tissue viability is $>$ 40%	No Category

211

212

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 9 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

213 **2. MATERIALS**

214 **2.1 LabCyte CORNEA-MODEL24**

215 **2.1.1 LabCyte CORNEA-MODEL24 KIT COMPONENTS**

216 LabCyte CORNEA-MODEL24 kit components are shown in **Table 2**.

217 **Table 2** LabCyte CORNEA-MODEL24 Kit Components

Component	Qty	Description
LabCyte CORNEA-MODEL24 plate	1 plate	Contains 24 culture inserts with tissues fixed in nutritive agar medium for transport (usable area: 0.3 cm ²).
Assay Medium	1 bottle	Basic medium for incubation (30 mL). Keep refrigerated.
24-well plate	1 plate	Blank plate for use in assay. Store at room temperature.

218


219 **2.1.2 SHIPMENT OF LabCyte CORNEA-MODEL24**

220 LabCyte CORNEA-MODEL24 is packed in a special Icompo container available from and delivered
221 by NIPPON EXPRESS CO., LTD. After the kit is delivered, examine the contents and make sure that all
222 components (LabCyte CORNEA-MODEL24 plate, assay medium, and 24-well assay plate) are included
223 in the package. Also confirm lot numbers and expiration dates. Record details in Methods Documentation
224 Sheet (MDS) 1. (See **MSD-1**).

225 NIPPON EXPRESS will pick up the Icompo container at a later date (generally, the day after delivery),
226 and it should be returned together with an shipping invoice and the insulating materials.

227

228

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 10 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

229 **2.1.3 INSTRUCTIONS FOR USE OF LabCyte CORNEA-MODEL24**

230 Always incubate all of the culture inserts after opening the package. Do not store the culture inserts
231 after opening.

232 The human corneal epithelial tissue cells used in the LabCyte CORNEA-MODEL24 originate from a
233 normal human donor and are HIV-, HBV-, HCV-, and HPV-negative. They are to be handled,
234 nevertheless, with due care and in accordance with laboratory biosafety guidelines for handling
235 human-derived materials.

236

237 **2.2 CONSUMABLES**

238 The following consumables are required.

239 *The following quantities are necessary to assay between one and six 6 test chemicals at a time.

- 240 • Assay Medium, 100 mL (J-TEC: 402250) 1
- 241 • Cell Counting Kit-8, 500 test (Dojindo: CK04) 4
- 242 • 24-well assay plate (Becton, Dickinson and Company: 353047) 7-8
- 243 • 96-well plate (Becton, Dickinson and Company: 353072 or equivalents) 1
- 244 • PBS, 500 mL (Invitrogen: 14190-144 or equivalents) 2 or 3 bottle
- 245 • Earle balanced salt solution (EBSS), 500 mL (SIGMA-ALDRICH: E3024) 1 bottle
- 246 • Sterile cotton buds (JAPAN COTTON BUDS: 10A754D or equivalents) 1 box
- 247 • Micro-pipette tips (sterile: 10~200 μ L, 200-1000 μ L)
- 248 • Microtubes (1.5mL)
- 249 • Dish (10cm)
- 250 • Paper towel

251


252 Convenient consumable items are shown followings.

253 Also, it would be convenient to have the following.

- 254 • Capillary & piston for positive-displacement-type pipette (10-100 μ L)

255

256

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 11 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

257 **2.3 OTHERS**

258 **2.3.1 EQUIPMENT/INSTRUMENTS**


- 259 • Safety cabinet (or clean bench)
- 260 • Water bath (37°C)
- 261 • CO₂ incubator (37°C, 5% CO₂, capable of maintaining high humidity)
- 262 • Autoclave
- 263 • 96-well multi-plate reader (required filters: 450 nm, 650 nm)
- 264 • Precision balance (0.1 mg)
- 265 • Aspirator
- 266 • Stop-watches
- 267 • Adjustable micro-pipette (10–200 µL, 200–1000 µL)
- 268 • Sharp-edged forceps (sterile)
- 269 • Micro spatula (sterile)
- 270 • Beaker (1–2 L: sterile)
- 271 • Sterilizable poly wash bottle (500–1000 mL: sterile) with wide mouth (mouth > 3-mm dia.)
- 272 • Mortar with pestle

273

274 Also, it would be convenient to have the following.

- 275 • Positive-displacement-type pipette (10–100 µL)

276

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 12 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

277 **3. TEST METHOD**

278 *Procedures described in Sections 3.1.1 to 3.1.3 and Sections 3.3.1 to 3.3.3 are to be performed
279 aseptically in a safety cabinet or clean bench. Procedures other than those mentions in the previous
280 sentence need not be performed aseptically. Refer to **Section 2.1.3 “INSTRUCTIONS FOR USE**
281 **OF LabCyte CORNEA-MODEL24”**.

282
283 **3.1 PREPARATIONS**

284 **3.1.1 POSITIVE CONTROL**

- 285 (1) Ethanol is used as a positive control for liquid test chemicals.
286 (2) Lauric acid is used as a positive control for solid test chemicals.


287
288 **3.1.2 NEGATIVE CONTROL**

- 289 (1) PBS is used as a negative control for liquid test chemicals.
290 (2) Non-treatment is a negative control for solid test chemicals.

291
292 **3.1.3 POLY WASH BOTTLE FOR PBS**

- 293 (1) Sterilize poly wash bottle by autoclave.
294 (2) Fill the sterilized poly wash bottle with sterile PBS.

295
296
297

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 13 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

298 **3.2 TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT**

299 There are two kinds of test chemicals that interfere with the WST-8 assay.

300 (a) Test chemicals that stain corneal epithelial tissues directly.

301 (b) Test chemicals that react directly with WST-8.

302 The test chemicals that can dye the corneal epithelial tissues will be extracted from the colored tissue
303 during WST-8 reaction, affecting the OD measurement. Also, the test chemicals that can directly
304 reduce the WST-8 medium will affect the OD measurement due to unexpected reduction reaction by
305 residues of test chemicals in the culture inserts. These test chemicals are needed to perform
306 additional experiments. The procedure is described below.

307

308

309 **3.2.1 DETECTION OF THE CHEMICALS THAT STAIN THE TISSUE**

310 **3.2.1.1 STEP 1 (PRELIMINARY TEST)**

311 (1) Add 50 μ L (Liquid) or 10mg (Solid) of the test chemical into wells of 24-well assay plate
312 preliminarily filled with 0.5mL of distilled water. Untreated distilled water is used as control.

313 (2) Close the lid of 24-well assay plate and incubate the mixture in CO₂ incubator for 15 minutes.

314 (3) After incubation, mix gently and evaluate the color change of the distilled water by visual check.

315 (4) When the color of the solution changes significantly, the test chemical is presumed to have the
316 potential to stain the tissue. A functional check on viable tissues (Step2) should be performed. If the
317 color of the solution does not change significantly, the test chemical is determined not to have a
318 potential to stain the tissue.


319 (5) Record the details of steps 1 to 4 above in **MDS 1-2**.

320

321 **3.2.1.2 STEP 2 (FUNCTIONAL CHECK ON VIABLE TISSUE)**

322 (1) Add 50 μ L (liquid) or 10 mg (solid) of the test chemical, which could clearly change the color of the
323 distilled water in step 1, onto the surface of the epidermis tissues. PBS are used as negative control.

324 (2) Go to the Section 3.3 EXECUTION OF THE TEST and perform the experiments according to the
325 procedures, expect for the WST-8 reaction. In the section of WST-8 assay, use the EBSS that does

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 14 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

326 not contain WST stock solution, instead of diluted WST-8 medium to evaluate the affected OD
327 values by colored diluted water. (correcting tissue).

328 (3) The corrected OD is calculated using the following formula.

329 **Corrected OD = A – (B – C)**

330

331 A : the OD of viable tissue exposed to a test chemical.

332 B : the mean OD of correcting tissue exposed to a test chemical.

333 C : the mean OD of correcting tissue exposed to the negative control.

334

335 (4) If a corrected OD is below 0, the OD is considered to be 0.

336 (5) When a cell viability that is calculated according to the procedures described in **Section 3.3.4.4** is
337 $\leq 40\%$, the test chemical is predicted to be an irritant (GHS Category 1 or 2) and there is no need
338 to calculate a corrected value.

339

340

341 **3.2.2 DETECTION OF CHEMICALS THAT DIRECTLY REDUCE WST-8**

342 **3.2.2.1 STEP 3 (PRELIMINARY TEST)**

343 (1) Dilute the cell counting kit-8 (WST-8 stock solution) with EBSS (Cell Counting Kit-8:EBSS = 1:10),
344 and then prepare the diluted WST-8 medium.


345 *Prepare WST-8 medium before use.


346 Dispense 0.3 mL of diluted WST-8 medium into each well of the 24-well assay plate.

347 (2) Add 50 μ L of a liquid test chemical or 10 mg of a test chemical to the wells of 24-well assay plate.
348 The diluted WST-8 medium is used as control.

349 (3) Put on the lid of 24-well assay plate and incubate in a CO₂ incubator for about 4 hours.

350 (4) After incubation, shake the mixture gently and evaluate the color change of the diluted WST-8
351 medium by visual check.

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 15 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.



Ethanol 5% Triton X100 Glycolic acid m-Phenylene diamine Monoethanolamine Calcium thioglycollate

352 **Photo 1** Example of test chemicals to interfere the WST-8 assay directly (STEP 1).

353 When test chemicals like m-Phenylene diamine, Monoethanolamine, or Calcium thioglycollate have
354 colored the diluted WST-8 medium, Step 4 must be performed.

355 (5) Significant coloring of the diluted WST-8 medium by the test chemical indicates the interference
356 between test chemicals and the WST-8 assay medium. The additional functional check is required for
357 these interference chemicals. Go to step 4 described below.

358 (6) Record the details of steps 1 to 5 above in **MDS 1-3**.

359

360 **3.2.2.2 STEP 4 (FUNCTIONAL CHECK ON FREEZE-KILLED TISSUE)**

361 (1) Add 50 µL of a liquid test chemical or 10 mg of solid a test chemical that clearly changed the color
362 of the diluted WST-8 medium (**3.2.2.1. STEP 3**) to the surface of the corneal epithelial tissues.

363 (2) Go to **Section 3.3 “EXECUTION OF THE TEST”** and perform the experiments using
364 freeze-killed corneal epithelial tissues, instead of using viable corneal epithelial tissues. . The
365 freeze-killed tissues are prepared by placing untreated Labcyte CRNEA-MODEL in -80°C or lower
366 for more than 30 min.

367 Record the details of freeze-killing the tissue in **MDS 1-3**.

368 (3) The corrected OD is calculated using the following formula.

369 **Corrected OD = A – (B – C)**

370


371 A : the OD of viable tissue exposed to a test chemical.

372 B : the mean OD of freeze-killed tissue (correcting tissue) exposed to a test chemical.

373 C : the mean OD of freeze-killed tissue (correcting tissue) exposed to the negative control.

374


375 (4) If a corrected OD is below 0, the OD is considered to be 0.

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 16 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

376 (5) Calculate a cell viability according to the procedures described in **Section 3.3.4.4**. In case that the
377 cell viability is $\leq 40\%$, the test chemical is predicted to be an irritant (GHS Category 1 or 2) and
378 there is no need to calculate a corrected value.

379

380

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 17 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

381 **3.3 EXECUTION OF THE TEST**

382

383 **3.3.1 PREPARATION OF LabCyte CORNEA-MODEL24 (DAY -1)**

384 (1) Warm the assay medium to 37°C for 30 minutes
385 in a water bath.

386 (2) Dispense 0.5 mL/well of the warm assay
387 medium into the six wells of the 1st row of each
388 assay plate for LIQUID/SOLID .

389 **See Fig. 1.**

390 (3) Open the LabCyte CORNEA-MODEL24
391 aluminum package.

392 (4) Open the LabCyte CORNEA-MODEL24 plate lid and pick up the culture inserts using sterile
393 forceps.

394 *Do not touch the surface of the corneal epithelial tissue in the culture inserts.

395 *Use forceps to remove any agar medium
396 sticking to the outside of the culture inserts.

397 (5) Transfer the culture inserts to the assay medium
398 in the wells of the 1st row using sterile forceps.

399 **See Fig. 2.**

400 (6) * Avoid the formation of air bubbles under the
401 tissue inserts. Close the lid on the plate and
402 place it in a CO₂ incubator.



403 (7) Incubate overnight (15–30 hours) until ready to perform **Section 3.3.2 “APPLICATION OF**
404 **LIQUID TEST CHEMICALS, RINSING AND POST-INCUBATION”** or **Section 3.3.3**
405 **“APPLICATION OF SOLID TEST CHEMICALS AND RINSING”**.


406 (8) Record the details of steps 1 to 7 above in MDS 2.

407

408

409



Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 18 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

410 **3.3.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1)**


411

412 **3.3.2.1 PREPARATION OF WELLS FOR POST-INCUBATION (3RD ROW)**

413 (1) Warm the assay medium to 37°C for 30 minutes
414 using a water bath.

415 (2) Take out the assay plate for LIQUID from the
416 CO₂ incubator.

417 (3) Open the lid of the assay plate for LIQUID, and
418 use a micropipette to fill the six wells of the 3rd
419 row with 0.5 mL/well of the warm assay medium.



420 **See Fig. 3.**


421 (4) Close the lid of the assay plate and perform **Section 3.3.2.2 “APPLICATION OF LIQUID TEST**
422 **CHEMICALS AND RINSING”** continuously.

423 (5) If the test chemicals are not applied immediately, store the Assay Plate for LIQUID in a CO₂
424 incubator until ready apply but for no more than 12 hours.

425 (6) Record the details of steps 1 to 5 above in **MDS 3-1**.

426


427

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 19 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

428 **3.3.2.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING**

429 (1) Take out the assay plate for LIQUID from the CO₂ incubator.

430 (2) Using a micropipette, apply 50 µL of a liquid test
431 chemical to the surface of the corneal epithelial
432 tissues in the 1st row of the assay plate. Each test
433 chemical is to be tested in three wells (N=3).
434 Carefully apply the test chemical to the central
435 part of each corneal epithelial tissue. After
436 application, close the lid of the assay plate and tap
437 the sides of the plate to spread the liquid test
438 chemicals to spread out over the entire corneal




439 epithelial surface. If necessary, use a micro spatula to spread the liquid test chemical over the entire
440 surface. Take care not to press down on the surface of the corneal epithelial with the spatula.

441 *For viscous LIQUID test chemicals, use a wide orifice cell saver tip (See **Photo 2.**) or positive
442 displacement type pipette.

443 *Use a pipette or other equipment to familiarize
444 yourself beforehand with the characteristics
445 the test chemicals.

446 *Assay no more than two test chemicals on one
447 24-well assay plate.




448 **See Fig. 4.**

449 Each chemical is tested in three wells, using
450 three tissues (N = 3).

451 (3) Apply test chemicals to each well at an interval of one to three minutes.

452 (4) Close the lid and incubate each well for 60 ±10 seconds in the cabinet at room temperature.

453 *Keep the lid of the assay plate closed at all times except when applying test chemicals. Leaving the
454 lid open could affect the quantity of the test chemical in the well due to air circulation in the
455 cabinet.

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 20 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

456 (5) Wait for 60 ± 10 seconds after applying a test chemical,
457 then open the assay plate for LIQUID and pick up a
458 culture insert with sterile forceps.

459 (6) Discard the test chemical on the tissue by decantation
460 and tapping it on a beaker. Fill the culture insert to
461 overflowing with PBS from a poly wash bottle.

462 * Keep PBS flush during rinsing to wash away the test
463 chemical from the tissue surface..

464 **See Photo 3.**

465 *To avoid damaging the tissue with too forceful a
466 stream, use a wide-mouth nozzle on the poly wash
467 bottle.

468 (7) Decant the PBS into the beaker. Remove as much of the
469 PBS inside the culture insert as possible by tapping it on
470 the beaker.

471 **See Photo 4.**

472 (8) Repeat steps 6 and 7 at least 10 times to remove as much as possible of the residual test chemical on
473 the tissue surface.

474 *Depending upon the physical properties of the test
475 chemical, bubbles might form in an insert during
476 washing. Continue washing until all bubbles
477 disappear.

478 (9) Using a sterile cotton bud, gently remove as much as
479 possible of the leftover PBS both inside and outside the
480 culture insert.

481 **See Photo 5.**

482 *Take care not to press down on the surface of the tissue
483 with the sterile cotton bud.

Photo 3 Rinse 1



Photo 4 Rinse 2






Photo 5 Rinse 3



Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 21 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

484

485 (10) Wipe the culture insert and then place it in the well
486 of the same column in the 3rd row.

487 **See Fig. 5.**

488 *Take care to prevent air bubbles from forming
489 under the culture inserts.

490 (11) Repeat steps 1 to 11 for all culture inserts at one- to
491 three-minute intervals.

492 (12) Record the details of steps 1 to 12 above in **MDS 3-1**.

493

494 3.3.2.3 POST-EXPOSURE INCUBATION

495 (1) After performing **Section 3.3.2.2 “APPLICATION OF LIQUID TEST CHEMICALS AND**
496 **RINSING”**, close the lid of the assay plate for LIQUID and place it in a CO₂ incubator as soon as
497 possible.


498 (2) Incubate for 24 ±1 hours.


499 (3) Record the details of steps 1 and 2 above in **MDS 3-1**.

500

501

502



Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 22 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

503 **3.3.3 APPLICATION OF SOLID TEST CHEMICALS AND RINSING (DAY 0–1)**

504


505 **3.3.3.1 PREPARATION OF WELLS FOR SOLID TEST CHEMICAL APPLICATION (2ND**
506 **ROW) AND FOR POST-INCUBATION (3RD ROW)**

507 (1) Warm the assay medium to 37°C for 30 minutes using a water bath.

508 (2) Take out the assay plate for SOLID from the CO₂
509 incubator.

510 (3) Open the lid of the assay plate for SOLID and use
511 a micropipette to fill the 12 wells in the 2nd and 3rd
512 rows with 0.5 mL/well of warm assay medium .

513 **See Fig. 6.**



514 (4) Close the lid of the assay plate and go to **Section**

515 **3.3.3.2 “APPLICATION OF SOLID TEST CHEMICALS”** continuously.

516 (5) If the test chemicals are not applied immediately, store the Assay Plate for LIQUID in a CO₂
517 incubator until ready apply but for no more than 12 hours.

518 (6) Record the details of steps 1 to 5 above in **MDS 3-2**.

519

520 **3.3.3.2 APPLICATION OF SOLID TEST CHEMICALS**

521 (1) Take out the assay plate for SOLID from the CO₂
522 incubator.

523 Using a precision balance, weigh out 10 ±2 mg of
524 the solid test chemicals. If necessary, crush and
525 grind the solid test chemicals in a mortar with
526 pestle. Apply the solid test chemical to the surface
527 of the corneal epithelial tissue.

528 **See Photo 6.**



529 If necessary, use a micro spatula to spread the test chemical gently over the entire surface. Use three
530 wells per test chemical (N = 3).

531 See Fig. 7.

532

533 (2) Place the exposed culture insert in the well of
534 the same column in the 2nd row (chemical
535 application).

536 See Fig. 8.

537 *Take care to prevent air bubbles from forming
538 under the culture inserts.

539 (3) Close the lid of the assay plate for SOLID and
540 place it in a CO₂ incubator. Incubate for
541 24 ±1 hours.

542 (4) Record the details of steps 1 to 3 above in
543 **MDS 3-2**.

544

545 3.3.3.3 RINSING OF SOLID TEST CHEMICALS

546 (1) After incubation, take out the assay plate for SOLID from the CO₂ incubator.

547 (2) Open the assay plate for SOLID and pick up a culture insert with sterile forceps.



548 (3) Discard test chemicals on the tissue by tilting the insert and tapping if on a beaker. Fill the culture
549 insert to overflowing with PBS from a poly wash bottle.

550 * Keep PBS flush during rinsing to wash away the test
551 chemical from the tissue surface.

552 See Photo 7.

553 *To avoid damaging the tissue with too forceful a
554 stream, use a wide-mouth nozzle on the poly wash
555 bottle.

556 (4) Tilt the insert to discard the PBS into the beaker.
557 Remove as much of the PBS inside the culture insert as
558 possible by tapping it on the beaker.



544

545 3.3.3.3 RINSING OF SOLID TEST CHEMICALS

546 (1) After incubation, take out the assay plate for SOLID from the CO₂ incubator.

547 (2) Open the assay plate for SOLID and pick up a culture insert with sterile forceps.

548 (3) Discard test chemicals on the tissue by tilting the insert and tapping if on a beaker. Fill the culture
549 insert to overflowing with PBS from a poly wash bottle.

550 * Keep PBS flush during rinsing to wash away the test
551 chemical from the tissue surface.

552 See Photo 7.

553 *To avoid damaging the tissue with too forceful a
554 stream, use a wide-mouth nozzle on the poly wash
555 bottle.

556 (4) Tilt the insert to discard the PBS into the beaker.
557 Remove as much of the PBS inside the culture insert as
558 possible by tapping it on the beaker.



559 **See Photo 8.**
560 (5) Repeat steps 3 and 4 at least 10 times to remove as
561 much as possible of the residual test chemical on the
562 tissue surface.

563 (6) Using a sterile cotton bud, gently remove as much as
564 possible of the leftover PBS both inside and outside
565 the culture insert.

566 **See Photo 9.**
567 (7) If it proves difficult to remove completely all the
568 residual test chemical from the corneal epithelial tissue
569 surface, remove as much as possible and continue to
570 step 8.

571 (8) Place the rinsed culture insert in the well of the same
572 column in the 3rd row (rinse).



573 **See Fig. 9.**

574 *Take care to prevent air bubbles from forming
575 under the culture inserts.

576 (9) Record the details of steps 1 to 8 above in **MDS**
577 **3-2.**

578 After rinsing step, perform **Section 3.3.4**
579 **“WST-8 assay”** continuously.

580
581



582 **3.3.4 WST-8 ASSAY (DAY 1)**

583

584 **3.3.4.1 PREPARATION OF WELLS FOR WST-8 ASSAY**

- 585 (1) Warm the EBSS to 37°C for 30 minutes using a water bath.
- 586 (2) Dilute the cell counting kit-8 (WST-8 stock solution) with EBSS (Cell Counting Kit-8:EBSS = 1:10),
- 587 and then warm the diluted WST-8 medium. Prepare the additional dilute WST-8 medium for blanks
- 588 in the WST-8 assay.

589 *Prepare WST-8 medium before performing the

590 WST-8 assay.

- 591 (3) Take out the assay plate for LIQUID from the CO₂
- 592 incubator or prepare the assay plate for SOLID.

- 593 (4) Open the lid of the assay plate and use a
- 594 micropipette to fill each well of the 4th row with
- 595 0.3 mL/well of the warm diluted WST-8 medium.

596 See Fig. 10.

597 Close the lid of the assay plate and perform **Section 3.3.4.2 “WST-8 REACTION”** continuously.



- 598 (5) Record the details of steps 1 to 4 above in **MDS 4-1**.

599

600 **3.3.4.2 WST-8 REACTION**

- 601 (1) Add about 20 mL of PBS each to two dishes: PBS dish 1 and PBS dish 2.
- 602 (2) Open an assay plate (for either LIQUID or SOLID) and pick up a culture insert with sterile forceps.
- 603 (3) Remove residual culture medium on a culture
- 604 insert by washing in PBS dish 1 (first time) and the
- 605 in PBS dish 2 (second time). After washing, wipe
- 606 the bottom with a paper towel.
- 607 (4) After washing and wiping, place the culture insert
- 608 in the well of the same column in the 4th row.

609 See Fig. 11.



- 610 * Avoid the formation of air bubbles under the culture inserts.
- 611 (4) Close the lid of the assay plate and place it in the CO₂ incubator.
- 612 Incubate for 4 hours ± 20 minutes.
- 613 Record the details of steps 1 to 5 above in **MDS 4-1**.

614

615 **3.3.4.3 SAMPLING THE REACTED WST-8 MEDIUM**

- 616 (1) After incubation, take out the assay plate from the CO₂ incubator.
- 617 (2) Open the lid of the assay plate and remove the culture inserts from the 4th row with forceps.
- 618 Transfer 200 µL of the reacted WST-8 dilution medium into the wells of a 96-well plate.
- 619 *Figs. 12A and 12B show typical allocations on a 96-well plate for both living and correcting tissue
- 620 (freeze-killed tissue or tissue reacted with only EBSS that does not contain WST-8 stock solution).

Fig.12A Allocation on a 96-well plate for living tissue

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank											
B	NC for LIQUID-1	NC for LIQUID-2	NC for LIQUID-3	PC for LIQUID-1	PC for LIQUID-2	PC for LIQUID-3	NC for SOLID-1	NC for SOLID-2	NC for SOLID-3	PC for SOLID-1	PC for SOLID-2	PC for SOLID-3
C	LIQUID 1-1	LIQUID 1-2	LIQUID 1-3	LIQUID 2-1	LIQUID 2-2	LIQUID 2-3	SOLID 1-1	SOLID 1-2	SOLID 1-3	SOLID 2-1	SOLID 2-2	SOLID 2-3
D	LIQUID 3-1	LIQUID 3-2	LIQUID 3-3	LIQUID 4-1	LIQUID 4-2	LIQUID 4-3	SOLID 3-1	SOLID 3-2	SOLID 3-3	SOLID 4-1	SOLID 4-2	SOLID 4-3
E	LIQUID 5-1	LIQUID 5-2	LIQUID 5-3	LIQUID 6-1	LIQUID 6-2	LIQUID 6-3	SOLID 5-1	SOLID 5-2	SOLID 5-3	SOLID 6-1	SOLID 6-2	SOLID 6-3
F	LIQUID 7-1	LIQUID 7-2	LIQUID 7-3	LIQUID 8-1	LIQUID 8-2	LIQUID 8-3	SOLID 7-1	SOLID 7-2	SOLID 7-3	SOLID 8-1	SOLID 8-2	SOLID 8-3
G	LIQUID 9-1	LIQUID 9-2	LIQUID 9-3	LIQUID 10-1	LIQUID 10-2	LIQUID 10-3	SOLID 9-1	SOLID 9-2	SOLID 9-3	SOLID 10-1	SOLID 10-2	SOLID 10-3
H												

621

622

Fig.12B Allocation on a 96-well plate for correcting tissue.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank											
B	NC for LIQUID-1	NC for LIQUID-2	NC for LIQUID-3				NC for SOLID-1	NC for SOLID-2	NC for SOLID-3			
C	LIQUID 1-1	LIQUID 1-2	LIQUID 1-3	LIQUID 2-1	LIQUID 2-2	LIQUID 2-3	SOLID 1-1	SOLID 1-2	SOLID 1-3	SOLID 2-1	SOLID 2-2	SOLID 2-3
D	LIQUID 3-1	LIQUID 3-2	LIQUID 3-3	LIQUID 4-1	LIQUID 4-2	LIQUID 4-3	SOLID 3-1	SOLID 3-2	SOLID 3-3	SOLID 4-1	SOLID 4-2	SOLID 4-3
E	LIQUID 5-1	LIQUID 5-2	LIQUID 5-3	LIQUID 6-1	LIQUID 6-2	LIQUID 6-3	SOLID 5-1	SOLID 5-2	SOLID 5-3	SOLID 6-1	SOLID 6-2	SOLID 6-3
F	LIQUID 7-1	LIQUID 7-2	LIQUID 7-3	LIQUID 8-1	LIQUID 8-2	LIQUID 8-3	SOLID 7-1	SOLID 7-2	SOLID 7-3	SOLID 8-1	SOLID 8-2	SOLID 8-3
G	LIQUID 9-1	LIQUID 9-2	LIQUID 9-3	LIQUID 10-1	LIQUID 10-2	LIQUID 10-3	SOLID 9-1	SOLID 9-2	SOLID 9-3	SOLID 10-1	SOLID 10-2	SOLID 10-3
H												

Add reacted WST-8 medium to the same locations as living tissues (Fig.12A).

623
624
625
626
627
628
629
630
631
632
633
634
635
636

- (3) Record the details of steps 1 to 2 above in **MDS 4-2**.


3.3.4.4 OPTICAL DENSITY MEASUREMENTS OF THE REACTION MEDIUM

- (1) Using a 96-well plate reader, measure OD at 450 nm and 650 nm and then used the following equation to determine a composite OD for each well.

$$\text{Composite OD} = (\text{OD}_{\text{TC}} \text{ at } 450 \text{ nm} - \text{OD}_{\text{blank}} \text{ at } 450 \text{ nm}) - (\text{OD}_{\text{TC}} \text{ at } 650 \text{ nm} - \text{OD}_{\text{blank}} \text{ at } 650 \text{ nm})$$

*If the plate reader can be programmed to perform this calculation automatically, then only the composite OD value need be recorded.

- (2) Calculate the mean OD for the negative control, a cell viability for each individual tissue, and a mean cell viability (including SD) for each test chemical using the following equations.

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 28 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

$$\text{Mean OD}_{\text{NC}} = \frac{\text{Sum of the OD}_{\text{NC}} \text{ for three replicate tissues}}{3}$$

637

$$\text{Tissue cell viability (\%)} = \frac{\text{Each tissue OD}_{\text{TC}}}{\text{Mean OD}_{\text{NC}}} \times 100$$

638


$$\text{Mean cell viability (\%)} = \frac{\text{Sum total of cell viability (\%)} \text{ for three replicate tissues}}{3}$$

639

640 (3) Record the details of steps 1 and 2 above in **MDS 4-2**.

641

642

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 29 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

643 **4. ASSESSMENT**

644

645 **4.1 CONDITIONS FOR A VALID TEST**

646 An eye irritation test is considered valid if all three of the following criteria have been met.

647

- 648 • Tissue viability: $0.5 \leq \text{mean OD (A450/650) measured value for negative control} \leq 1.6$
- 649 • Positive control: mean tissue viability for positive control $\leq 40\%$
- 650 • SD: SD (negative control, positive control and each test chemicals) of tissue viability of 3
- 651 identically treated replicates $\leq 18\%$

652

653 **4.2 ASSAY CRITERIA**

654 The criteria for in vitro prediction are shown below.

655 After exposure to a chemical, if cell viability is 40% or less, the chemical is predicted to be an
656 irritant (GHS Category 1 or 2), otherwise it is predicted to be a non-irritant (GHS No Category).

657 See Table 3.


658

659 **Table 3** Prediction model of LabCyte24 EIT

Tissue Viability	Prediction
$\leq 40\%$	Category 1 or 2
$> 40\%$	No Category

660

661


Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 30 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

662 **5. References**

663 Green H. (1978) Cyclic AMP in relation to proliferation of the epidermal cell: New view. Cell, 15,
664 801-811.

665 Rheinwald J.G. and Green H. (1975) Serial cultivation of strains of human epidermal keratinocytes:
666 The formation of keratinizing colonies from single cells. Cell, 6, 331-343.

667

Version 2.7		Page 31 of 45
October, 2018	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	 Japan Tissue Engineering Co., Ltd.

668 **MDS 1-1:**
669 **RECEIPT OF LabCyte CORNEA-MODEL24 (2.1.2)**

670
671 Laboratory name: _____ Test name: _____ Test no. : _____

672
673 1. LabCyte CORNEA-MODEL24
674

675 Date received: _____
676

677 Lot no.: _____
678

679 Expiration date: _____
680

681
682 Accessories: Assay medium, 30mL (Lot no.: _____ Expiration date: _____)
683

684 24 well assay plate
685
686

<u>Note</u>

687
688
689 2. Assay medium
690

691
692 Date received: _____
693


694 Lot no.: _____
695

696 Expiration date: _____
697

<u>Note</u>

699
700
701 Date: _____ Operator: _____ Check date: _____ Study director: _____
702

703
704 Secretariat Check date: _____ Name: _____
705
706

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 32 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

708
709
710
711
712
713
714
715
716
717
718
719
720
721
722
723
724
725
726
727
728
729
730

MDS 1-2:
TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT
STEP-1 (3.2.1.1)

Laboratory name:_____ Test name: _____ Test no. : _____

1. Add distilled water (0.5 mL) to the wells of the 24-well assay plate.
To add distilled water (0.5 mL) Execution date/time: _____
2. Apply test chemicals to the wells of the 24-well assay plate.
3. Culture the 24-well assay plate in CO₂ incubator for 4 hours.
Time of WST-8 reaction started: _____
Time of WST-8 reaction completed: _____
4. Check the color of water.
5. Test chemical information and check list of coloring potential.

Test chemical	Physical state	Amount	Coloring	Test chemical.	Physical state	Amount	Coloring
PBS (NC)	LIQUID			Non treatment (NC)			
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	

731
732

Note

733
734
735
736
737
738
739
740

Date:_____ Operator:_____ Check date:_____ Study director:_____

Secretariat Check date:_____ Name:_____

741 **MDS 1-3:**
742 **TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT**
743 **STEP-3 (3.2.2.1)**

744 Laboratory name: _____ Test name: _____ Test no. : _____

746 1. Preparation of WST-8 dilution medium

747 Warm EBSS for 30 minutes. Time/date: _____

748 CCK-8: (Lot no.: _____ Expiration date: _____)

749 EBSS: (Lot no.: _____ Expiration date: _____)

750 Volume _____ mL Time/date completed: _____

751 2. Add WST-8 dilution medium (0.3mL) to the wells of the 24-well assay plate.

752 To add WST-8 dilution medium (0.3mL) Time/date executed: _____

753 3. Apply test chemicals to the wells of the 24-well assay plate.

754 4. Culture the 24-well assay plate in CO₂ incubator for 4 hours.

755 Time of WST-8 reaction started: _____

756 Time of WST-8 reaction completed: _____

757 5. Check the color of WST-8 medium.


758 6. Test chemical information and checked list of WST-8 assay interfere.

Test chemical	Physical state	Amount	WST-8 assay interfere	Test chemical.	Physical state	Amount	WST-8 assay interfere
PBS (NC)	LIQUID			Non treatment (NC)			
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	

773 Note

774 Date: _____ Operator: _____ Check date: _____ Study director: _____

775 Secretariat Check date: _____ Name: _____

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 34 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

783
784
785
786
787
788
789
790
791
792
793
794
795
796
797
798
799
800
801
802
803
804
805
806
807
808
809
810
811
812
813
814
815
816

MDS 1-4:
PREPARATION OF FREEZE KILLED TISSUE (3.2.2)

Laboratory name: _____ Test name: _____ Test no. : _____

1. Transfer LabCyte CORNEA-MODEL24 tissues to 50 mL tube or appropriate sterile container.
2. Freeze tissues in the -80°C deep-freezer for 30 minutes (1st freezing).
Store for 30 minutes. Time/date: _____
3. Thaw tissues in the 37°C incubator for 15 minutes.
Store for 15 minutes. Time/date: _____
4. Freeze tissues in the -80°C deep-freezer for more than 30 minutes (2nd freezing).
Store for more 30minutes. Time/date: _____
5. Just before using, thaw tissues in the 37°C incubator for 15 minutes.
Store for 15 minutes. Time/date: _____

Note

817
818
819
820
821
822
823
824

Date: _____ Operator: _____ Check date: _____ Study director: _____

Secretariat Check date: _____ Name: _____

863 **MDS 3-1(LIQUID):**

864 **APPLICATION OF LIQUID TEST CHEMICALS, RINSING AND POST-INCUBATION (3.3.2)**

865 Laboratory name: _____ Test name: _____ Test no. : _____

- 866
867 1. Warm up the assay medium and add 0.5mL of the assay medium to the wells of the 3rd row on the
868 24-well assay plate for LIQUID.

869 Assay medium: (Lot no.: _____ Expiration date: _____)

870 Warm for 30 minutes. Time/date: _____

871 Add 0.5mL of assay medium to each well Time/date: _____ Number of plate: _____

- 872
873 2. Apply test chemicals to the LabCyte CORNEA-MODEL24.

874 Time/date execution started: _____

- 875 3. LIQUID test chemical information

Test chemical.	Lot no.	Physical state	Test chemical vol.	Time of application	Time of rinsing	Exposure period (1 minute)
PBS (Negative control)		LIQUID	50 µL	1	:	<input type="checkbox"/>
				2	:	<input type="checkbox"/>
				3	:	<input type="checkbox"/>
Ethanol (Positive control)		LIQUID	50 µL	1	:	<input type="checkbox"/>
				2	:	<input type="checkbox"/>
				3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>

- 880
881 4. After exposure to test chemical for 1 minute, wash out the LabCyte CORNEA-MODEL24 and
882 transfer the culture inserts to the 3rd row on the 24-well assay plate.

883 PBS: (Lot no.: _____ Expiration date: _____)


884 Hit PBS stream on the tissue surface directly.

885 Confirm that there are no bubbles under the cell culture insert. Time/date completed: _____

- 886
887 5. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 24 hours.

888 Time/date post-incubation started: _____

889 Time/date post-incubation completed: _____

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 37 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.


893

<u>Note</u>

894
895
896
897
898
899
900
901
902
903
904
905
906

Date: _____ Operator: _____ Check date: _____ Study director: _____

Secretariat Check date: _____ Name: _____

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 38 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

907
908
909
910
911
912
913
914
915
916
917
918
919
920
921
922
923
924
925
926
927
928
929
930
931
932
933
934
935
936
937
938
939
940
941
942

MDS 3-2(SOLID):
APPLICATION OF SOLID TEST CHEMICALS AND RINSING (3.3.3)

Laboratory name: _____ Test name: _____ Test no. _____

1. Warm up the assay medium and add 0.5mL of the assay medium to the wells of 2nd and 3rd row on the 24-well assay plate for SOLID.

Assay medium: (Lot no.: _____ Expiration date: _____)

Warm for 30 minutes. Time/date: _____
Add 0.5mL of assay medium to each well Time/date: _____ Number of plate: _____

2. Apply test chemicals to the LabCyte CORNEA-MODEL24 and transfer the culture inserts to the 2nd row on the 24-well assay plate.

Time/date started: _____ Time/date completed: _____

3. SOLID test chemical information

Test chemical	Lot no.	Physical state	Crush and grind	Test chemical amount.	Time of application	Exposure period (24hours)
Non treatment (Negative control)					:	<input type="checkbox"/>
Lauric acid (Positive control)		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>


4. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 24 hours.

Time/date exposure started: _____
Time/date exposure completed: _____

5. After exposure to test chemical, wash out the LabCyte CORNEA-MODEL24 and transfer the culture inserts to the 3rd row on the 24-well assay plate.

PBS: (Lot no.: _____ Expiration date: _____)
Hit PBS stream on the tissue surface directly.

Time/date started: _____ Time/date completed: _____


Version 2.7		Page 39 of 45
October, 2018	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	 Japan Tissue Engineering Co., Ltd.

Note

943
944
945
946
947
948

Date: _____ Operator: _____ Check date: _____ Study director: _____

Secretariat Check date: _____ Name: _____

Version 2.7	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 40 of 45
October, 2018		 Japan Tissue Engineering Co., Ltd.

949 **MDS 4-1:**
950 **WST-8 ASSAY (3.3.4.1, 3.3.4.2)**

951
952 Laboratory name: _____ Test name: _____ Test no. : _____

953
954 1. Preparation of WST-8 dilution medium

955 Warm EBSS for 30 minutes. Time/date: _____
956
957 CCK-8: (Lot no.: _____ Expiration date: _____)
958 EBSS: (Lot no.: _____ Expiration date: _____)
959 Volume _____ mL Time/date completed: _____
960

961
962 2. Add WST-8 dilution medium (0.3mL) to the wells in the 4th row on the 24-well assay plate.

963 To add WST-8 dilution medium (0.3mL) Time/date executed: _____
964
965

966
967 3. After the operation, the blotted tissue transfer to wells of 4th row of 24-well assay plate.

968 Time/date started: _____ Time/date completed: _____
969
970

971 Confirm that there are no bubbles under the cell culture insert.

972
973
974 4. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 4 hours.

975 Time of WST-8 reaction started: _____
976
977 Time of WST-8 reaction completed: _____
978
979

980
981
982
983
984
985
986
987

<u>Note</u>

988
989
990
991
992 Date: _____ Operator: _____ Check date: _____ Study director: _____
993
994

995 Secretariat Check date: _____ Name: _____
996
997

999 **MDS 4-2:**
1000 SAMPLING OF REACTING WST-8 SOLUTION (3.3.4.3) AND MEASUREMENT (3.3.4.4)

1001
1002 Laboratory name: _____ Test name: _____ Test no. : _____
1003

1004 1. Reacting WST-8 solution (200 µL) is transferred to each well on the 96-well plate.

1005
1006 Transfer to the 96-well plate.
1007

1008 Time/date executed: _____
1009
1010

1011 Sample location on 96-well plate.

	<u>LIQUID</u>						<u>SOLID</u>					
	1	2	3	4	5	6	7	8	9	10	11	12
A	blank											
B	PBS-1	PBS-2	PBS-3	Ethanol-1	Ethanol-2	Ethanol-3	Non-treatment-1	Non-treatment-2	Non-treatment-3	Lauric acid-1	Lauric acid-2	Lauric acid-3
C												
D												
E												
F												
G												
H												

1012
1013
1014
1015 2. Analyze extract OD at 450nm and 650nm, and calculate the OD(450nm-650nm).
1016

1017 Analyze OD at 450nm and 650nm.

1018 Calculate the OD (450nm-650nm).

1019 Calculate cell viability and SD.

1020 Cell viability and SD are recorded on a separate data sheet.

1021 The data sheet is attached to the back of this sheet.


1022 Check for input errors.

1023 Time/date executed: _____
1024
1025

Note

1026
1027 Date: _____ Operator: _____ Check date: _____ Study director: _____
1028
1029


1030 Secretariat Check date: _____ Name: _____
1031

Version 2.7	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 42 of 45
October, 2018		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>


1032 **REVISION HISTORY**

1033


Rev.	CONTENT	Date Revised
Ver.1.1	1) First version	Aug., 2012
Ver.1.2	1) Revised clerical error.	Aug., 2012
Ver.1.3	1) Revised clerical error.	Sep., 2012
Ver.1.3PS	<p>Added the following supplementary explanation.</p> <p>1) Added more detail explanation about the conditions of WST-8 reaction.</p> <p>2) Added the formula of SD.</p>	Sep., 2012
Ver.2.1	<p>1) Revised clerical error.</p> <p>2) In the section 1.2 “BACKGROUND”, changed statement about animal testing (Draize test).</p> <p>3) In the section 3.3. “TEST METHOD”, explained the washing protocol of the <u>LIQUID</u> and <u>SOLID</u> test chemical more briefly.</p> <p>4) In the section 3.3. “TEST METHOD”, changed the WST-8 dilution rate with PBS from 1:10 to 1:5.</p> <p>5) In the section 3.3. “TEST METHOD”, changed the reaction period of WST-8 from 5 hours to 4 hours.</p> <p>6) In the section 3.3. “TEST METHOD”, changed the condition of WST-8 reaction from shaking to standing.</p> <p>7) In the section 3.3. “TEST METHOD”, changed the application-amount of <u>SOLID</u> chemicals from 50mg to 10mg.</p> <p>8) As assay acceptance criteria, added that SD (test chemicals) of tissue viability of 3 identically treated replicates $\leq 20\%$.</p> <p>9) In the prediction model of this EIT, changed the cut-off value of the mean viability from 50% to 40%.</p>	May, 2013

Version 2.7	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 43 of 45
October, 2018		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

Ver.2.2	<ol style="list-style-type: none"> 1) Revised clerical error. 2) In the section 3.3. "TEST METHOD", additionally explained that the solid test chemical is crush and grind in a mortar with pestle if necessary. 3) In the section 3.3.2 "APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1)", explained temperature condition (room temperature) at chemical application. 	Sep., 2013
Ver.2.3	<ol style="list-style-type: none"> 1) Revised clerical error. 2) In the section 3.3. "TEST METHOD", explained about a blank preparation at WST-8 assay. 3) In the section 3.3. "TEST METHOD", changed the condition of WST-8 reaction from standing to shaking. 5) In the section 3.3. "TEST METHOD", explained about a blank preparation at WST-8 assay. 	Feb., 2014
Ver.2.3.1	<ol style="list-style-type: none"> 1) Revised clerical error. 2) About the prediction result of eye irritation, changed the classification from irritation/no irritation to GHS classification. 	Mar., 2014
Ver.2.3.2MTT	<ol style="list-style-type: none"> 1) At the analysis of cell viability, changed the assay method from WST-8 assay to MTT assay. 	Jul., 2014
Ver.2.4.1	<ol style="list-style-type: none"> 1) At the analysis of cell viability, changed the assay method from MTT assay to WST-8 assay. 2) In the section 3.3. "TEST METHOD", changed the condition of WST-8 reaction from shaking to standing. 3) In the section 3.3. "TEST METHOD", changed the dilution solution of WST-8 reaction from PBS to EBSS. 4) In the section 3.3. "TEST METHOD", changed the reaction period of WST-8 from 5 hours to 4 hours. 5) The standard of the additional testing was mentioned about a 	Jan., 2015

Version 2.7	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 44 of 45
October, 2018		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

	<p>borderline result.</p> <p>9) In the prediction model of this EIT, changed the cut-off value of the mean viability from 40% to 50%.</p>	
Ver.2.4.2	1) Changed the assay acceptance criterion of negative control from 0.5 < and <2.0 to 0.6 < and <1.5.	Mar., 2015
Ver.2.4.2m	1) The judgement of the chemical which is a result of the borderline was added.	Jul., 2015
Ver.2.5.1	<p>1) Revised clerical error.</p> <p>2) In the prediction model of this EIT, changed the cut-off value of the mean viability from 50% to 40%.</p> <p>3) Changed the judgement of equivocal results.</p> <p>4) Changed the assay acceptance criterion of negative control from 0.6 < and <1.5 to 0.5 < and <1.3.</p> <p>5) In the section 3.2. "TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT", revised description has clearly.</p> <p>6) Added MDS1-2 and MDS 1-3.</p>	Aug., 2015
Ver.2.5.1m	<p>1) Revised clerical error.</p> <p>2) Changed allocation for a 96-well plate of pattern which consists of freeze killed tissue (Fig.12B).</p>	Sep., 2015
Ver. 2.5.1mr	1) Revised Fig.2B	Sep.,2015
Ver. 2.5.2	<p>1) In the section 3.3. "TEST METHOD", explained the washing protocol of the <u>LIQUID</u> and <u>SOLID</u> test chemical more briefly.</p> <p>2) In the MDS 3-1 and the MDS 3-2, added the check box about the washing procedure.</p> <p>3) In the MDS 3-2, added the check box about crush and grind of</p>	Sep.,2015

Version 2.7	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 45 of 45
October, 2018		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

	test chemicals.	
Ver. 2.5.4	<p>1) In the section 3.3. "TEST METHOD", added attention point of the washing procedure.</p> <p>2) In the section 4.3 "ASSESSMENT FLOWCHART", Change numbers of test run from three independent run to single run.</p>	June, 2016
Ver. 2.5.5	<p>1) Revised clerical error.</p> <p>2) As assay acceptance criteria, changed that SD (negative control, positive control and test chemicals) of tissue viability of 3 identically treated replicates from $\leq 20\%$ to $\leq 18\%$.</p>	October, 2016
Ver.2.5.6	<p>1) Revised clerical error.</p> <p>2) Detection protocol of coloring interference is changed from using WST-8 medium to distilled water and correction of coloring interference is changed from using freeze-killed tissue to using living tissue without WST-8 reaction.</p>	February, 2017
Ver.2.6	<p>1) Changed the assay acceptance criterion of negative control from $0.5 \leq$ and ≤ 1.3 to $0.5 \leq$ and ≤ 1.6.</p>	June, 2018
Ver.2.7	<p>1) Revised clerical errors. There is no change of protocol.</p>	

1034

1035