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In vitro Cytotoxicity Test Methods for Estimating Starting Doses
for Acute Oral Systemic Toxicity Test

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ADME	Absorption, distribution, metabolism, and elimination
ANOVA	Analysis of variance
ATC	Acute Toxic Class method
ATWG	Acute Toxicity Working Group
BRD	Background Review Document
	本報告書では <i>In Vitro</i> Cytotoxicity Test Methods for Estimating Acute Oral Systemic Toxicity (NIH Publication No. 07-4518)を示す。
CNS	Central nervous system
CV	Coefficient of variation
DMEM	Dulbecco's Modification of Eagle's Medium
DMSO	Dimethyl sulfoxide
D-PBS	Dulbecco's phosphate buffered saline
EC ₅₀	Concentration of a substance that produces 50% of the maximum possible response for that substance
ECBC	U.S. Army Edgewood Chemical Biological Center
ECVAM	European Centre for the Validation of Alternative Methods
ETOH	Ethanol (Ethyl alcohol)
FAL	FRAME Alternatives Laboratory
FDP	Fixed Dose Procedure
FRAME	Fund for the Replacement of Animals in Medical Experiments
GHS	Globally Harmonized System (of Classification and Labeling of Chemicals)

GLP	Good Laboratory Practices
HBSS	Hanks' balanced salt solution
IC ₅₀	Concentration producing 50% inhibition of the endpoint measured
ICCVAM	Interagency Coordinating Committee for the Validation of Alternative Methods
IIVS	Institute for In Vitro Sciences
i.p.	Intraperitoneal
i.v.	Intravenous
KBM	Keratinocyte basal medium
Kow	Octanol-water partition coefficient
LD ₅₀	Dose that produces lethality in 50% of test animals
MTT	3-(4, 5-Dimethylthiazol-2-yl)-2, 5-diphenyl tetrazolium bromide
N	Number (of substances)
NHK	Normal human epidermal keratinocytes
NICEATM	National Toxicology Program Center for the Evaluation of Alternative Toxicological Methods
NR	Neutral red
NRU	Neutral red uptake
OD	Optical density
OD ₅₄₀	Optical density (absorbance) at a wavelength of 540 nm
PBS	Phosphate buffered saline
PC	Positive control
pH	Power of hydrogen

QA	Quality assurance
QC	Quality control
r	Pearson correlation coefficient
R ²	Coefficient of determination
r _s	Spearman correlation coefficient
RC	Registry of Cytotoxicity
RI	Radioisotope
RTECS [®]	Registry of Toxic Effects of Chemical Substances
SD	Standard deviation
SLS	Sodium lauryl sulfate
SMT	Study management team
3T3	BALB/c mouse fibroblasts, clone A31 (ATCC # CCL-163)
UDP	Up-and-Down Procedure
VC	Vehicle control
ZEBET	Zentralstelle zur Erfassung und Bewertung von Ersatz- und Ergänzungsmethoden zum Tierversuch (German Center for Documentation and Evaluation of Alternative Methods to Animal Experiments)

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
RC millimole regression (BRD より抜粋)

Figure 1-1 RC Millimole Regression for *In Vitro* Cytotoxicity (IC₅₀) and Rat and Mouse Acute Oral LD₅₀ Values for 347 Chemicals

Abbreviations: RC=Registry of Cytotoxicity; IC₅₀=Geometric mean (of multiple endpoints and cell types) test substance concentration that reduces cell viability by 50%; LD₅₀=Dose producing death in 50% of the animals tested.

The heavy line shows the fit of the data to a linear regression model, $\log(LD_{50}) = 0.435 \times \log(IC_{50}) + 0.625$; $r=0.67$. The thinner lines show the empirical prediction interval ($\pm \log 5$, or ± 0.699) that is based on the anticipated precision for the prediction of LD₅₀ values from cytotoxicity data (Halle 1998, 2003).

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
423

OECD/OCDE

ANNEX 2c: TEST PROCEDURE WITH A STARTING DOSE OF 300 MG/KG BODY WEIGHT

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Phase Ia: Laboratory Evaluation

Development of a positive control database for each laboratory

Phase Ib: Laboratory Evaluation

Limited substance testing to demonstrate the reliability of the protocol

Phase II: Laboratory Qualification

Evaluation of protocol refinements

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Table 1-2 Regulatory Classification Systems for Acute Oral Toxicity

Regulatory Agency (Authorizing Act)	Animals	Endpoint	Classification
EPA (FIFRA)	Use current EPA or OECD protocol	Death ¹	I - LD ₅₀ ≤ 50 mg/kg II - 50 < LD ₅₀ ≤ 500 mg/kg III - 500 < LD ₅₀ ≤ 5000 mg/kg IV - LD ₅₀ > 5000 mg/kg
CPSC (Federal Hazardous Substances Act)	White rats, 200-300 g	Death ¹ within 14 days for ≥ half of a group of ≥ 10 animals	Highly toxic - LD ₅₀ ≤ 50 mg/kg Toxic - 50 mg/kg < LD ₅₀ < 5 g/kg
OSHA (Occupational Safety and Health Act)	Albino rats, 200-300 g	Death ¹ , duration not specified.	Highly toxic - LD ₅₀ ≤ 50 mg/kg Toxic - 50 < LD ₅₀ < 500 mg/kg
DOT (Federal Hazardous Material Transportation Act)	Male and female young adult albino rats	Death ¹ within 14 days of half the animals tested. Number of animals tested must be sufficient for statistically valid results.	Packing Group I - LD ₅₀ ≤ 5 mg/kg Packing Group II - 5 < LD ₅₀ ≤ 50 mg/kg Packing Group III - LD ₅₀ < 500 mg/kg (liquid) LD ₅₀ < 200 mg/kg (solid)
OECD Guidance for Use of GHS (2001b)	Protocols not specified	Not specified	I - LD ₅₀ ≤ 5 mg/kg II - 5 < LD ₅₀ ≤ 50 mg/kg III - 50 < LD ₅₀ ≤ 300 mg/kg IV - 300 < LD ₅₀ ≤ 2000 mg/kg V - 2000 < LD ₅₀ ≤ 5000 mg/kg Unclassified - LD ₅₀ > 5000 mg/kg

Abbreviations: EPA=U.S. Environmental Protection Agency; OECD=Organisation for Economic Co-operation and Development; LD₅₀=Dose producing death in 50% of the animals tested; CPSC=U.S. Consumer Product Safety Commission; FIFRA=Federal Insecticide, Fungicide, and Rodenticide Act; OSHA=U.S. Occupational Safety and Health Administration; DOT=U.S. Department of Transportation; GHS=Globally Harmonized System of Classification and Labelling of Chemicals (UN 2005).

¹Guidance Document on the Recognition, Assessment and Use of Clinical Signs as Humane Endpoints for Experimental Animals Used in Safety Evaluation calls for humane killing of moribund animals (OECD 2000). Moribund animals that are humanely euthanized are accepted as deaths.

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Abbreviations: ECBC=Edgewood Chemical Biological Center; FAL=Fund for the Replacement of Animals in Medical Experiments Alternatives Laboratory; IIVS=Institute for *In Vitro* Sciences; 3T3=BALB/c 3T3 fibroblasts; N=Number of substances used to calculate the regression; NHK=Normal human epidermal keratinocytes; NRU=Neutral red uptake; R²=Coefficient of determination.

¹Log IC₅₀ in mM; log LD₅₀ in mmol/kg.

²Regression based on a single point per substance (i.e., the geometric mean of the within laboratory replicate IC₅₀ values and the reference rat acute oral LD₅₀ from Table 4-2).

³Regression based on a single point per substance (i.e., the geometric mean of the geometric mean IC₅₀ values obtained for each laboratory and the reference rat acute oral LD₅₀ from Table 4-2).

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Table 6-3 Outlier Substances for the RC and the 3T3 and NHK NRU Methods When the RC Millimole Regression is Used¹

Substances Included in the RC Identified as Outliers in:		
RC ²	3T3 ³	NHK ⁴
	Acetaminophen (+)	
	Arsenic III trioxide (-)	Arsenic III trioxide (-)
		Aminopterin (-)
5-Aminosalicylic acid (+)		5-Aminosalicylic acid (+)
Busulfan (-)	Busulfan (-)	Busulfan (-)
Caffeine (-)		Caffeine (-)
Cycloheximide (-)	Cycloheximide (-)	Cycloheximide (-)
Dibutyl phthalate (+)	Dibutyl phthalate (+)	Dibutyl phthalate (+)
	Diethyl phthalate (+)	Diethyl phthalate (+)
Digoxin (-)	Digoxin (-)	
Disulfoton (-)	Disulfoton (-)	Disulfoton (-)
Epinephrine bitartrate (-)	Epinephrine bitartrate (-)	Epinephrine bitartrate (-)
Ethanol (+)	Ethanol (+)	Ethanol (+)
Lindane (-)	Lindane (-)	
Mercury II chloride (-)	Mercury II chloride (-)	Mercury II chloride (-)
		Methanol (+)

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Abbreviations: RC=Registry of Cytotoxicity; 3T3=BALB/c 3T3 fibroblasts; NHK=Normal human epidermal keratinocytes; NRU=Neutral red uptake; (-)=Toxicity was underpredicted by the IC₅₀ and RC millimole regression (i.e., the LD₅₀ value predicted by the IC₅₀ was higher than the *in vivo* LD₅₀ value); (+)=Toxicity was overpredicted by the IC₅₀ and RC millimole regression (i.e., the LD₅₀ value predicted by the IC₅₀ was lower than the *in vivo* rodent LD₅₀ value).

[Note: Empty cells indicate that the substance was not an outlier for that particular IC₅₀ value.]

¹Log LD₅₀ (mmol/kg) = 0.435 log IC₅₀ (mM) + 0.625. Log LD₅₀ (mmol/kg) values for outlier substances were >0.699 from the RC millimole regression.

²Using RC IC₅₀ in the RC millimole regression for the 58 RC substances tested in the validation study.

³Using the 3T3 NRU IC₅₀ in the RC millimole regression for the 70 reference substances that yielded IC₅₀ values from any laboratory in the validation study.

⁴Using the NHK NRU IC₅₀ in the RC millimole regression the RC for the 71 reference substances that yielded IC₅₀ values from any laboratory in the validation study.

Bolded substances have active metabolites *in vivo* (see Table 3-7).

Substances that showed evidence of insolubility (i.e., precipitates) during testing (see Table 5-11) are identified by italics.

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Table 6-7 Prediction of GHS Acute Oral Toxicity Category by the 3T3 and NHK NRU Test Methods and the RC Rat-Only Millimole Regression¹

Reference Rat Oral LD ₅₀ ^a (mg/kg)	3T3 -Predicted GHS Category (mg/kg)						Total	Accuracy	Toxicity Over-predicted	Toxicity Under-predicted
	LD ₅₀ < 5	5 < LD ₅₀ ≤ 50	50 < LD ₅₀ ≤ 300	300 < LD ₅₀ ≤ 2000	2000 < LD ₅₀ ≤ 5000	LD ₅₀ > 5000				
LD ₅₀ < 5	0	2	0	4	0	0	6 ^d	0%	0%	100%
5 < LD ₅₀ ≤ 50	0	1	6	3	1	0	11 ^e	9%	0%	91%
50 < LD ₅₀ ≤ 300	0	0	5	7	0	0	12	42%	0%	58%
300 < LD ₅₀ ≤ 2000	0	1	2	13	0	0	16	81%	19%	0%
2000 < LD ₅₀ ≤ 5000	0	0	0	10	0	0	10 ^f	0%	100%	0%
LD ₅₀ > 5000	0	0	0	8	2	2	12 ^g	17%	83%	0%
Total	0	4	13	45	3	2	67	31%	34%	34%
Predictivity	0%	25%	38%	29%	0%	100%				
Category Overpredicted	0%	25%	15%	40%	67%	0%				
Category Underpredicted	0%	50%	4%	31%	33%	0%				
Reference Rat Oral LD ₅₀ ^a	NHK -Predicted Toxicity Category (mg/kg)						Total	Accuracy	Toxicity Over-predicted	Toxicity Under-predicted
	LD ₅₀ < 5	5 < LD ₅₀ ≤ 50	50 < LD ₅₀ ≤ 300	300 < LD ₅₀ ≤ 2000	2000 < LD ₅₀ ≤ 5000	LD ₅₀ > 5000				
LD ₅₀ < 5	0	1	2	3	0	0	6 ^d	0%	0%	100%
5 < LD ₅₀ ≤ 50	0	2	5	3	1	0	11 ^e	18%	0%	82%
50 < LD ₅₀ ≤ 300	0	1	6	5	0	0	12	50%	8%	42%
300 < LD ₅₀ ≤ 2000	0	1	2	12	1	0	16	75%	19%	6%
2000 < LD ₅₀ ≤ 5000	0	0	0	10	0	0	10 ^f	0%	100%	0%
LD ₅₀ > 5000	0	0	0	7	6	0	13 ^g	0%	100%	0%
Total	0	5	15	40	8	0	68	29%	40%	31%
Predictivity	0%	40%	40%	30%	0%	0%				
Category Overpredicted	0%	40%	13%	43%	75%	0%				
Category Underpredicted	0%	20%	47%	28%	25%	0%				

Abbreviations: GHS=Global Harmonized System of Classification and Labelling of Chemicals (UN 2002); 3T3=BALB/c 3T3 fibroblasts; NHK=Normal human epidermal keratinocytes; NRU=Neural red uptake; RC=Registry of Cytotoxicity. Shaded cells are those containing the correct predictions.

¹The RC rat-only millimole regression is $\log LD_{50} \text{ (mmol/kg)} = \log IC_{50} \text{ (mM)} \times 0.439 + 0.621$. Numbers in table represent numbers of substances.

²Reference rat oral LD₅₀ values in mg/kg from Table 4-2.

³Epinephrine bitartrate excluded because no rat reference acute oral LD₅₀ was identified (see Table 4-2).

⁴Colchicine excluded because no rat acute oral LD₅₀ was identified (see Table 4-2).

⁵Carbon tetrachloride excluded because no laboratory attained sufficient toxicity for the calculation of an IC₅₀.

⁶Methanol excluded because no laboratory attained sufficient toxicity for the calculation of an IC₅₀.

⁷Propylparaben excluded because no rat acute oral LD₅₀ was identified (see Table 4-2).

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Table 10-3 Animal Use¹ for the UDP² by GHS Acute Oral Toxicity Category³ Using Starting Doses Based on the 3T3 and NHK NRU Test Methods with the RC Rat-Only Millimole Regression⁴

GHS Acute Oral Toxicity Category ³	Number of Reference Substances	Dose-mortality Slope = 2.0			Dose-mortality Slope = 8.3		
		With Default Starting Dose ⁵	With IC ₅₀ -Based Starting Dose ⁶	Animals Saved ⁷	With Default Starting Dose ⁵	With IC ₅₀ -Based Starting Dose ⁶	Animals Saved ⁷
3T3 NRU Test Method							
LD ₅₀ ≤ 5 mg/kg	6	11.32 ± 0.20	10.19 ± 0.70	1.14 (10.0%)	9.70 ± 0.28	8.74 ± 0.43	0.96 (9.9%)
5 < LD ₅₀ ≤ 50 mg/kg	11	9.68 ± 0.23	9.74 ± 0.45	-0.07 (-0.7%)	8.46 ± 0.28	8.54 ± 0.47	-0.08 (-1.0%)
50 < LD ₅₀ ≤ 300 mg/kg	12	7.76 ± 0.10	8.18 ± 0.21	-0.42 (-5.5%)	6.61 ± 0.19	6.90 ± 0.19	-0.29 (-4.3%)
300 < LD ₅₀ ≤ 2000 mg/kg	16	8.53 ± 0.21	8.14 ± 0.21	0.38 (4.5%)	7.46 ± 0.24	7.15 ± 0.19	0.31* (4.1%)
2000 < LD ₅₀ ≤ 5000 mg/kg	10	10.73 ± 0.10	9.46 ± 0.15	1.28* (11.9%)	9.17 ± 0.23	7.96 ± 0.31	1.21* (13.2%)
LD ₅₀ > 5000 mg/kg	12	9.87 ± 0.34	8.29 ± 0.49	1.58* (16.0%)	7.76 ± 0.59	6.18 ± 0.69	1.58* (20.3%)
NHK NRU Test Method							
LD ₅₀ ≤ 5 mg/kg	6	11.21 ± 0.24	10.47 ± 0.71	0.75 (6.7%)	9.66 ± 0.27	8.95 ± 0.52	0.71 (7.3%)
5 < LD ₅₀ ≤ 50 mg/kg	11	9.65 ± 0.16	9.99 ± 0.45	-0.34 (-3.5%)	8.43 ± 0.26	8.77 ± 0.49	-0.33 (-3.9%)
50 < LD ₅₀ ≤ 300 mg/kg	12	7.78 ± 0.11	8.12 ± 0.21	-0.34 (-4.4%)	6.57 ± 0.19	6.85 ± 0.19	-0.28 (-4.2%)
300 < LD ₅₀ ≤ 2000 mg/kg	16	8.55 ± 0.22	8.03 ± 0.23	0.52* (6.1%)	7.49 ± 0.25	7.00 ± 0.20	0.49* (6.5%)
2000 < LD ₅₀ ≤ 5000 mg/kg	10	10.75 ± 0.08	9.54 ± 0.20	1.21* (11.3%)	9.17 ± 0.23	8.06 ± 0.29	1.11* (12.1%)
LD ₅₀ > 5000 mg/kg	13	9.87 ± 0.32	8.41 ± 0.44	1.47* (14.8%)	7.66 ± 0.59	6.18 ± 0.69	1.47* (19.2%)

Abbreviations: 3T3=BALB/c 3T3 fibroblasts; GHS=Globally Harmonized System of Classification and Labelling of Chemicals (UN 2005); NHK=Normal human epidermal keratinocytes; RC=Registry of Cytotoxicity; UDP=Up-and-Down Procedure.

*Statistically significant (p<0.05) by a one-sided Wilcoxon signed rank test. Percentage difference shown in parentheses.

¹Mean numbers of animals used ± standard errors for 10,000 simulations for each substance with an upper limit dose of 5000 mg/kg. Although the simulations used whole animals, averaging the results over a large number of simulations produced fractional numbers. Results are provided for 67 substances in the 3T3 NRU test method and 68 substances in the NHK NRU test method. Substances were categorized using the rat acute oral LD₅₀ reference values in mg/kg from Table 4-2.

²OECD (2001a); EPA (2002a).

³UN (2005).

⁴The RC rat-only millimole regression is $\log LD_{50} (\text{mmol/kg}) = 0.439 \log IC_{50} (\text{mM}) + 0.621$.

⁵Default starting dose = 175 mg/kg.

⁶The starting dose was one default dose lower than the predicted LD₅₀ calculated using the IC₅₀ value for each reference substance in the RC rat-only millimole regression. The IC₅₀ value for each reference substance was randomly selected from the distribution of values obtained during the testing with each method.

⁷Difference between mean animal use with the default starting dose and mean animal use with the predicted starting dose.

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Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM), National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods (NICEATM), National Institute of Environmental Health Sciences, National Institutes of Health, U. S. Public Health Service, Department of Health and Human Services, BACKGROUND REVIEW DOCUMENT, *in vitro* Cytotoxicity Test Methods for Estimating Acute Oral Systemic Toxicity, NIH Publication No: 07-4518

Borenfreund E, Puerner JA. 1985. Toxicity determination *in vitro* by morphological alterations and neutral red absorption. *Toxicol Lett* 24:119-124.

Clemedson C, Barile FA, Ekwall Ba, Gon. *Toxicol Lett* 24:119-124. K, et al. 1998a. MEIC evaluation of acute systemic toxicity. Part III. *In vitro* results from 16 additional methods used to test the first 30 reference chemicals and a comparative cytotoxicity analysis. *Altern Lab Anim* 26 (suppl 1):93-129.

Clemedson C, Andersson M, Aoki Y, Barile FA, Bassi AM, Calleja MC, et al. 1998b. MEIC evaluation of acute systemic toxicity. Part IV. *In vitro* results from 67 toxicity assays used to test reference chemicals 31-50 and a comparative cytotoxicity analysis. *Altern Lab Anim* 26 (suppl 1):131-183.

Ekwall B. 1983. Screening of toxic compounds in mammalian cell cultures. *Ann New York Acad Sci* 407:64-77.

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Halle W. 2003. The Registry of Cytotoxicity: Toxicity testing in cell cultures to predict acute toxicity (LD₅₀) and to reduce testing in animals. *Altern Lab Anim* 31:89-198. (English translation of Halle 1998)

ICCVAM. 2001a. Report Of The International Workshop On *In vitro* Methods For Assessing Acute Systemic Toxicity. NIH Publication No. 01-4499. Research Triangle Park, NC:National Institute

for Environmental Health Sciences. Available: <http://iccvam.niehs.nih.gov/> [accessed 01 November 2006].

ICCVAM. 2001b. Guidance Document On Using In vitro Data To Estimate In vivo Starting Doses For Acute Toxicity. NIH Publication No. 01-4500. Research Triangle Park, NC: National Institute for Environmental Health Sciences. Available:<http://iccvam.niehs.nih.gov/> [accessed 01 November 2006].

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OECD. Guideline for Testing of Chemicals, 425, Acute Oral Toxicity – Up-and-Down Procedure.

OECD. Guideline for Testing of Chemicals, 420, Acute Oral Toxicity – Fixed Dose Method.

OECD. Guideline for Testing of Chemicals, 423, Acute Oral Toxicity – Acute Toxic Class Method.

OECD Environmental Health and Safety Publications Series on Testing and Assessment No. 19
Guidance Document on The Recognition, Assessment, and Use of Clinical Signs as Humane Endpoints for Experimental Animals Used in Safety Evaluation